



Bringing Nature and Technology Together

Standard Reference Solution @ 25°C

for calibration of
pH Meters

Specimen

INSTRUCTIONS

Pour an adequate quantity of reference solution into a clean container. Immerse the pH electrode in pH reference solution for 1 minute before calibrating. Rinse the electrode with distilled water after each pH reading to avoid contamination of the reference solution. Never pour used solution back into this bottle.

pH 4.01

**KEEP OUT OF REACH
OF CHILDREN**

Manufactured by
 **GH**
GENERAL HYDROPONICS
PO Box 1576, Sebastopol, CA 95473, USA
genhydro.com



946 ml (32 fl. oz)

s072506am